


+ Introducing ALLEVYN[®] Ag+ Antimicrobial Foam Dressing



Smith+Nephew

ALLEVYN[®] Ag+
Antimicrobial Foam Dressing


www.smith-nephew.com/allevyn

Indications

ALLEVYN Ag+ Dressing is indicated for the management of exuding wounds such as leg and foot ulcers, pressure ulcers, traumatic and surgical wounds, superficial and partial thickness burns. ALLEVYN Ag+ Dressing can also be used under compression.



**LOW - HIGH
ABSORBENCY**



Features and benefits



Superior[†] antimicrobial efficacy

ALLEVYN[®] Ag+ Antimicrobial Foam Dressing is superior[†] to Mepilex[™] Border Ag in reducing bacterial bioburden in the dressing.⁶



Helps protect the skin

ExuLOCK Advanced Lock-in Technology to minimize leakage of fluid.¹



Helps support daily patient activities

ALLEVYN Ag+ Dressing is showerproof.⁷



HighFLEX Technology

Flexible and conformable to optimize retention and patient comfort during movement.^{8,9}




ComfortSTAY Technology

Silicone adhesive is suitable for fragile skin and helps minimize pain and trauma on removal.⁵



ALLEVYN Ag+ Dressing has a 7-day wear time (5 days on sacrum)



ALLEVYN Ag+ Antimicrobial Foam Dressing

S+N Code	Size	Box
66027662	ALLEVYN Ag+ 3in x 3in	10
66027663	ALLEVYN Ag+ 4in x 4in	10
66027664	ALLEVYN Ag+ 5in x 5in	10
66027665	ALLEVYN Ag+ 7in x 7in	10
66027669	ALLEVYN Ag+ Sacrum 6.8in x 6.9in	10
66027670	ALLEVYN Ag+ Sacrum 8.5in x 9.1in	10

*Tested *in vitro*. [†]Tested *in vivo* using n=1 batch of each dressing. At day 1; p=0.006, day 3; p=0.003, day 7; p=0.002.

Advanced Wound Management
Smith & Nephew, Inc
Fort Worth, TX 76109
USA

www.smith-nephew.com
T 800-876-1261
F 727-392-6914

[®]Trademark of Smith+Nephew
All Trademarks acknowledged
©2025 Smith+Nephew
47841_US

For detailed product information, including indications for use, contraindications, warnings and precautions, please consult the product's Instructions for Use (IFU) prior to use.

References

1. Smith+Nephew. Internal report CSD.AWM.24.024. 2. Smith+Nephew. Internal Report. CSD.AWM.24.026. 3. Smith+Nephew. Internal Report. CSD.AWM.24.031. 4. Smith+Nephew. Internal report CSD.AWM.24.051. 5. Smith+Nephew. Internal report EA/AWM/GENERAL/001v2. 6. Smith+Nephew. Internal Report CSD.AWM.24.014. 7. Smith+Nephew. Internal report CSD. AWM.24.027. 8. Smith+Nephew. Internal report CSD.AWM.24.025. 9. Smith+Nephew. Internal report CSD.AWM.24.030.